

## Datasheet

### Ibuprofen polyclonal antibody

**Catalog Number:** PAB14735

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against Ibuprofen.

**Immunogen:** Ibuprofen conjugated with carrier protein.

**Host:** Rabbit

**Applications:** IHC

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Lyophilized

**Purification:** Ammonium sulfate precipitation

**Recommend Usage:** Immunohistochemistry

(1:1000-1:5000)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from PBS, pH 7.4

**Storage Instruction:** Store at 4°C on dry atmosphere.

After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years.

Aliquot to avoid repeated freezing and thawing.