

Datasheet

Pseudomonas aeruginosa **polyclonal antibody**

Catalog Number: PAB14743

Regulation Status: For research use only (RUO)

Product Description: Rat polyclonal antibody raised against *Pseudomonas aeruginosa*.

Immunogen: *Pseudomonas aeruginosa*

Host: Rat

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Recommend Usage: ELISA (1:2000-1:10000)

Western Blot (1:1000-1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from PBS (50% glycerol)

Storage Instruction: Store at 4°C on dry atmosphere.

After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years.

Aliquot to avoid repeated freezing and thawing.