

Datasheet

EIF2B3 mouse monoclonal antibody (hybridoma)

Catalog Number: H00008891-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant EIF2B3.

Immunogen: EIF2B3 (NP_065098.1, 1 a.a. ~ 452 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MEFQAVVMAVGGGSRMTDLTSSIPKLLPVGNKPLIW
YPLNLLERVGFEEVIVVTTTRDVQKALCAEFKMKMKPDI
VCIPDDADMGTADSLRYIYPKLKTDVLVLSCDLITDVAL
HEVVDFLFRAYDASLAMLMRKGQDSIEVPVGQKGKKKA
VEQRDFIGVDSTGKRLLFMANEADLDEELVIKGSILQK
HPRIRFHTGLVDAHLYCLKKYIVDFLMENGSISSIRSELI
PYLVRKQFSSASSQQGQEEKEEDLKKKELKSLDIYSFI
KEANTLNLAPYDACWNACRGDRWEDLSRSQVRCYVH
IMKEGLCSRSTLGLYMEANRQVPKLLSALCPEEPPV
HSSAQIVSKHLVGVDSLIGPETQIGEKSSIKRSVIGSSC
LIKDRVTITNCLLMNSVTVEEGSNIQGSVICNAVIEKG
ADIKDCLIGSGQRIEAKAKRVNEVIVGNDQLMEI

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Entrez GeneID: 8891

Gene Symbol: EIF2B3

Gene Alias: EIF-2B, EIF2Bgamma