

Datasheet

ZMYM3 mouse monoclonal antibody (hybridoma)

Catalog Number: H00009203-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant ZMYM3.

Immunogen: ZMYM3 (AAH13009.1, 1 a.a. ~ 495 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MDPSDFPSPFDPLTLPEKPLAGDLPVDMFEGEDLLES
QTAPTRGWAPPGPSPSSGALDLLDTPAGLEKDPGVLD
GATELLGLGGLLYKAPSPPEVDHGPETLAWDAGDQ
TLEPGPGGQTPEVVPPDPGAGANSCSPEGLLEPLAPD
SPITLQSPHIEEEETTSIATARRGSPGQEEELPQQGPQ
SPNAPPSPSVGETLGDGINSSQTKPGGSSPPAHPSLP
GDGLTAKASEKPPERKRSEVRRAEPPKPEVVDSTES
IPVSEDESDAMVDDPNDEDFVPRPRRSPRMSLRSSV
SQRAGRSVAGTKMTCAHCRTPLQKGQTAYQRKGLPQ
LFCSSSCLTTFSKKPSGKKTCTFCKKEIWNTKDSVVAQ
TGSGGSFHEFCTSVCLSLYEAQQQRPIQSGDPADAT
RCSICQKTGEVLHEVSNNGSVVHRLCSDSCFSKFRANK
GLKTNCCDQCGAYIYTKTGSPGPELLFHEGQQKRFCN
TTCLGAYKKVGPPE

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Entrez GeneID: 9203

Gene Symbol: ZMYM3

Gene Alias: DXS6673E, KIAA0385, MYM, XFIM,
ZNF198L2, ZNF261