

## Datasheet

### BZW1 mouse monoclonal antibody (hybridoma)

**Catalog Number:** H00009689-M

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody  
raised against a full-length recombinant BZW1.

**Immunogen:** BZW1 (NP\_055485.2, 1 a.a. ~ 353 a.a)  
full-length recombinant protein with GST tag. MW of the  
GST tag alone is 26 KDa.

**Sequence:**

MNNQKQQKPTLSGQRFKTRKRDEKERFDPTQFQDCII  
QGLTETGTDLEAVAKFLDASGAKLDYRRYAETLFDILV  
AGGMLAPGGTLADDMRTDVCVFAAQEDLETMQAFA  
QVFNKLIRRYKYLEKGFEDVKKLLFLKGFSESERNK  
LAMLTGVLLANGTLNASILNSLYNENLVKEGVSAFAV  
KLFKSWINEKDINAVAASLRKVSMDNRLMELFPANKQS  
VEHFTKYFTEAGLKELSEYVRNQQTIGARKELQKELQE  
QMSRGDPFKDIILYVKEEMKKNNIHFMAFQKIVVLFYK  
AEVLSEEPILKWYKDAHVAKGKSVFLEQMKKFVEWLK  
NAEEEESEAEEGD

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, WB-Re, WB-Tr  
(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Entrez GeneID:** 9689

**Gene Symbol:** BZW1

**Gene Alias:** BZAP45, KIAA0005, Nbla10236