

Datasheet

ACAP1 mouse monoclonal antibody (hybridoma)

Gene Symbol: ACAP1

Gene Alias: CENTB1, KIAA0050

Catalog Number: H00009744-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant ACAP1.

Immunogen: ACAP1 (NP_055531.1, 1 a.a. ~ 740 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MTVKLDFEECLKDSPRFRASIELVEAEVSELETRLEKLL
KLGTGLLESGRHYLAASRAFFVVGICDLARLGPPEPMM
AECLEKFTVSLNHLKLDSHAELLDATQHTLQQQIQTLVK
EGLRGFREARRDFWRGAESLEAALTHNAEVPRRRRAQ
EAEEAGAAALRTARAGYRGRALDYALQINVIEDKRKFDI
MEFVLRLVEAQATHFQQGHEELSRLSQYRKELGAQLH
QLVLNSAREKRDMEQRHVLLKQKELGGEEPEPSLRE
GPGGLVMEGHFLKRASNAFKTWSRRWFTIQSNQLVY
QKKYKDPVTVVVDDLRLCTVKLCPDSERRFCFEVVST
SKSCLLQADSERLLQLWVSAVQSSIASAFSQARLDDS
PRGPGQGSGHLAIGSAATLGSGGMARGREPGGVGH
VVAQVQSVDGNAQCCDCREPAPEWASINLGVTLCIQC
SGIHRSLGVHFSKVRSLTDSWEPELVKLMCELGNVII
NQIYEARVEAMAVKKPGPSCSRQEKEAWIHAKYVEKK
FLTKLPEIRGRRGGRPRGQPPVPPKPSIRPRPGSL
RSKPEPPSEDLGSLHPGALLFRASGHPPSLPTMADAL
AHGADVNWVNGGQDNATPLIQATAANSLACEFLLQN
GANVNQADSAGRGPLHHATILGHTGLACFLKRGADL
GARDSEGRDPLTIAMETANADIVTLLRLAKMREAAEQ
GQAGDETYLDIFRDFSLMASDDPEKLSRRSHDLHTL

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Entrez GeneID: 9744