

Datasheet

ZBTB39 mouse monoclonal antibody (hybridoma)

Gene Symbol: ZBTB39

Gene Alias: KIAA0352

Catalog Number: H00009880-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant ZBTB39.

Immunogen: ZBTB39 (NP_055645.1, 1 a.a. ~ 712 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MGMRIKLQSTNHPNNLLKELNKCRLSETMCDVTIVVG
SRSFPAHKAVLACAAGYFQNLFLNTGLDAARTYVVDI
TPANFEKVLVSFVYTSelfTDLINVGVIYEVAERLGMEDL
LQACHSTFPDLESTARAKPLTSTSESHSGTLSCPSAEP
AHPLGELRGGDYLGADRNYPVLPDAGGSYKEEEKN
VASDANHSLHLPQPPPPPKTEDHDTAPFTSIPSM
TQPLLGTVSTGIQTSTSSCQPYKVQSNQDFSKNSFLT
DNAVDITTGTNSCLSNSEHSKDPGFGQMDLQLEDLG
DDDQFEDPAEDIGTTEEVIELSDDSEDELAFGENDNR
ENKAMPCQVCKKVLEPNILIRQHARDHVDLLTGNCK
VCETHFQDRNSRVTHVLSHIGIFLSCDMCETKFFTQ
WQLTLHRRDGIFENNIIVHPNDPLPGKLGLFSGAASPE
LKCAACGKVLAKDFHVVRGHILDHLNLKGQACSVCDQ
RHLNLCSLMWHTLSHLGISVFSCSVCANSFVDWHLLE
KHMAVHQSLDALFHCRCSQSFKSEAARYHYVSQH
KCNGLDARPGFGLQHPALQKRKLPAEEFLGEELALQ
GQPGNSKYSCVKCGKRFAHTSEFNHYHRIHTGEKPY
QCKVCHKFFRGRSTIKCHLKTSGALMYRCTVCGHYS
STLNLMSKHVGVHKGSLPPDFTIEQTFMYIIHSKEADK
NPDS

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Entrez GeneID: 9880