

Datasheet

MYD88 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00004615-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human MYD88 protein.

Immunogen: MYD88 (NP_002459.1, 1 a.a. ~ 296 a.a) full-length human protein.

Sequence:

MAAGGPGAGSAAPVSSTSSLPLAALNMRVRRRLSLFL
NVRTQVAADWTALAEEMDFEYLEIRQLETQADPTGRL
LDAWQGRPGASVGRLLLELLTKLGRDDVLLELGPSIEE
DCQKYILKQQQEEAEKPLQVAADVSSVPRTAELAGITT
LDDPLGHMPERFDFICYCPSDIQFVQEMIRQLEQTN
YRLKLCVSDRDVLPGTCVWSIASSELIEKRCRRMVVVVS
DDYLSKECDFQTKFALSLSPGAHQKRLIPIKYKAMKK
EFPSILRFITVCDYTNPCTKSWFWTRLAKALSLP

Host: Rabbit

Reactivity: Human, Rat

Applications: IF, IP, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4615

Gene Symbol: MYD88

Gene Alias: MYD88D