

Datasheet

MYF6 purified MaxPab mouse polyclonal antibody (B02P)

Catalog Number: H00004618-B02P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human MYF6 protein.

Immunogen: MYF6 (NP_002460.1, 1 a.a. ~ 242 a.a) full-length human protein.

Sequence:

MMMDLFETGSYFFYLDGENVTLQPLEVAEGSPLYPGS
DGTLSQCQDQMPPEAGSDSSGEEHVLAPPGLQPPHC
PGQCLIWACKTCKRKSAPTDRRKAATLRERRRLKKIN
EAFEALKRRTVANPNQRLPKVEILRSAISYIERLQDLLH
RLDQQEKMQELGVDPFSYRPKQENLEGADFLRTCSS
QWPSVSDHSRGLVITAKEGGASIDSSASSSLRCLSSIV
DSISSEERKLPCVEEVVEK

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4618

Gene Symbol: MYF6

Gene Alias: MRF4, bHLHc4