



ZPACK13948098 PLASTIC BOTTLE LIQ 100 ML DIN 28



SPECIFICATIONS

Material :	HDPE
Neck Type :	DIN 28
Filling Opening :	28,00 mm
Bottom Diameter :	44,00 mm
Height :	96,00 mm

Remark : Packaging materials and specifications shown may be modified without preceding notice
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011
fax +32(0)16 385 385
info@be.vwr.com
www.vwr.com

Researchpark Haasrode 2020
Geldenaaksebaan 464 B-3001 Leuven
RPR/RPM Leuven
BTW/TVA BE 0403.593.343