

Datasheet

NO-Lysine polyclonal antibody

Catalog Number: PAB14738

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against NO-Lysine.

Immunogen: NO-Lysine conjugated with carrier protein.

Host: Rabbit

Applications: ELISA, ICC, IHC

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Purification: Ammonium sulfate precipitation

Recommend Usage: ELISA (1:2000-1:10000)

Immunohistochemistry (1:1000-1:5000)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from PBS, pH 7.4

Storage Instruction: Store at 4°C on dry atmosphere.

After reconstitution with 50 uL of distilled water and 50 uL of glycerol, store at -20°C for 2 years.

Aliquot to avoid repeated freezing and thawing.