

Datasheet

IL33 monoclonal antibody, clone 12B3C4

Catalog Number: MAB8086

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against recombinant IL33.

Clone Name: 12B3C4

Immunogen: Recombinant protein corresponding to
IL33.

Host: Mouse

Reactivity: Human, Mouse

Applications: ELISA, IHC, WB-Re
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Recommend Usage: Western Blot (1 ug/mL)
The optimal working dilution should be determined by
the end user.

Storage Buffer: In PBS (0.02% sodium azide)

Storage Instruction: Store at 4°C for three months. For
long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.